

Video OTOSCOPE

Let Our Image Enhance Your Image




Superior Image Quality
Is The Best Way To Describe
Our Video Oscope System

info@progresshearing.co.uk

Video Companion

Bringing Video Otoscopy To Your Computer



The remarkably clear resolution is particularly helpful in detecting excessive cerumen, foreign bodies, external otitis, PE tube status, cholesteatoma, TM perforations and performing cerumen management

- **Deluxe 3mm Otoscope Probe**
- **Compact Color Video Camera**
- **Focus Always Feature**
- **0-150 Watt Halogen Light Source**
- **USB 2.0 Video Capture Cable**
- **NOAH™ & TIMS® Compatible**

The battery operated LED light source is designed for the practitioner that requires portability



PROGRESS HEARING (UK)

info@progresshearing.co.uk

CE OS-S-ADH093-1