


# A Bold New Design With A Plus



Size:  
4.75" x 4.75" x 1.25"

Lighted Probe Mics  
Identify The Test Ear

**The AVANT A<sup>2</sup>D<sup>+</sup> & The  
AVANT REM Speech<sup>+</sup> Are  
The Complete Audiological  
Testing & Fitting System**

[www.progresshearing.co.uk](http://www.progresshearing.co.uk)

# AVANT The Ultimate Office<sup>+</sup>



Carrying Case Included

***A Complete Audiometric  
Testing & Fitting System  
with an Ergonomic Design,  
Integrated QuickSIN™ Testing  
plus DSL & NAL-NL1 Targets  
for Speech Mapping***

#0086   
OT-S-ADH111-1

**The AVANT A2D<sup>+</sup>**  
Is A Powerful Dual Channel  
Audiometer Featuring  
Integrated Quick-SIN  
Testing Capability

**The AVANT REM Speech<sup>+</sup>**  
Live Speech Mapping System  
Features the Sure-Probe™  
Lighted Microphone System  
Allowing For Monaural Or  
Simultaneous Binaural  
Measurements



**Our Sure-Probe™ offers an  
adjustable loop and probe  
retention design which allows for  
easy placement in any ear**

PROGRESS HEARING<sup>(UK)</sup>

**info@progresshearing.co.uk  
www.progresshearing.co.uk**